

Product Datasheet

RASSF3 antibody (orb19817)

Description	Goat polyclonal antibody to RASSF3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	RASSF3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 26kDa band observed in Human Heart lysates (calculated MW of 27.6kDa according to NP_835463.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	CKPPQTSPNSGK
Clonality	Polyclonal
MW	27.6
Entrez	283349
Dilution Range	ELISA: 1:2000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19817
(0.5ug/ml)
staining of
Heart ly...