

Product Datasheet

Sh2d1b / Sh2d1b2 (Mouse) antibody (orb19818)

Description

Goat polyclonal antibody to Sh2d1b /

Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Eat2A / Eat2B (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	ELNVYENTDEEYVD
Clonality	Polyclonal
MW	15.2
Dilution Range	ELISA: 1:1000
Expiration Date	12 months from date of receipt.