


Product Datasheet

CD56 Prediluted antibody (orb198233)

Description	Prediluted monoclonal antibody C5.9 in PBS containing 0.08% sodium azide
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Extract from KG1a cell line.
Concentration	See vial for concentration
Storage	Product should be stored at 4-8°C. DO NOT FREEZE
Note	For research use only
Application notes	Product provided in dropper bottle and should be used as provided. Use 2-3 drops per section or sufficient solution to cover section.
Isotype	IgG2b
Clonality	Monoclonal
Source	Derived from the hybridization of mouse Sp2/0 myeloma cells with spleen cells from BALB/c F1 mice immunized with KG1a cell line.
Uniprot ID	P13591
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pe...