


Product Datasheet

CD4, CD8 antibody (FITC, PE) [Out of stock] (orb198266)

Description	Mouse polyclonal antibody to CD4, CD8 (FITC, PE)
Species/Host	Gallus
Conjugation	HRP
Tested Applications	FACS
Preservatives	Phosphate Buffered Saline, 50% Glycerol with with 0.1% 2-chloroacetamide.
Concentration	See vial for concentration.
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles.
Note	For research use only
Application notes	Chickens were immunized with highly purified cervid IgG isolated from deer/elk serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbent and conjugated to horseradish peroxidase.
Isotype	IgG1 (F)/IgG1(PE)
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Peripheral
B...