

## Product Datasheet

### CD5, CD20 antibody (FITC, PE) [Out of stock] (orb198286)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Mouse Polyclonal antibody to CD5, CD20 (FITC, PE)  |
| <b>Species/Host</b>        | Gallus   |
| <b>Conjugation</b>         | HRP  |
| <b>Tested Applications</b> | FACS   |
| <b>Preservatives</b>       | Phosphate buffered saline, pH 7.2, 50% glycerol with 0.1% 2-chloroacetamide  |
| <b>Concentration</b>       | See vial for concentration.  |
| <b>Storage</b>             | Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles.  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Chickens were immunized with highly purified IgG from normal mouse serum. After multiple immunizations, eggs were collected and the IgY fraction prepared and the Fc specific antibodies immunoaffinity purified and conjugated to horseradish peroxidase. |
| <b>Isotype</b>             | IgG2a (F)/IgG1 (PE)  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |