

Product Datasheet

Lysophosphatidic Acid Acetyltransferase beta antibody (orb198343)

Description Rabbit Polyclonal antibody to Lysophosphatidic Acid Acetyltransferase beta

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications EIA, WB

Immunogen Synthetic peptide derived from the LPAAT-protein

Preservatives Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration See vial for concentration

Storage Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note For research use only


Application notes Antibody can be used for Western blotting (5-10 µg/ml) and EIA. Optimal concentration should be evaluated by serial dilutions.

Isotype IgG

Uniprot ID **O15120**

Hazard Information This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date 12 months from date of receipt.



Western blot analysis of human brain lys...