

Product Datasheet

TSGA2 antibody (orb19835)

Description	Goat polyclonal to TSGA2
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TSGA2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 37kDa band observed in Mouse Testis lysates (calculated MW of 35.1kDa according to Human NP_543136.1 and 34.2kDa according to Mouse NP_079566.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	HLNHRYQGKFLNKNP
Clonality	Polyclonal
MW	35.1
Entrez	89765
Dilution Range	ELISA: 1:4000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19835
(0.1ug/ml)
staining of
Mouse Te...