

## Product Datasheet

### Sphingosine Kinase 1 antibody (orb198363)

**Description** Provided as ligand affinity purified antibody in phosphate buffered saline with 0.08% sodium azide

**Species/Host** Rabbit

**Reactivity** Bovine, Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** ELISA, WB

**Concentration** See vial for concentration

**Storage** Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

**Note** For research use only

**Application notes** Antibody can be used for Western blotting (1:400 dilution) and ELISA. Other applications not yet tested. Optimal concentration should be evaluated by serial dilutions.

**Isotype** IgG

**Source** Rabbits were immunized with a synthetic peptide derived from the C-terminal of human Sphingosine Kinase 1

**Uniprot ID** [Q9NYA1](#)

**Hazard Information** This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

**Expiration Date** 12 months from date of receipt.



Western blot analysis of MDBK cell lysat...