


Product Datasheet

TXL2 antibody (orb19838)

Description	Goat polyclonal to TXL2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TXL2 / TXNDC6
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 38kDa band observed in Human Uterus lysates (calculated MW of 29.7kDa according to NP_835231.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Sadek et al, J Biol Chem. 2003 Apr 11;278(15):13133-42; 12569107). Recommended concentration: 1-3 µg/ml.
Protein Sequence	ESSTQPRLKITDLD
Clonality	Polyclonal
MW	29.7
Entrez	347736
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19838
(2ug/ml)
staining of
Human Uter...