

Product Datasheet

3T3 cells (orb198382)

Description	3T3 cells
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Form/Appearance	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Suggested use: A volume of 10-25ul is suggested for immunoblotting.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.