

## Product Datasheet

### Go 6983 (orb198422)

|                        |                                                                                                                                                                                                           |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Go 6983 is a pan-PKC inhibitor against for PKC $\alpha$ , PKC $\beta$ , PKC $\gamma$ and PKC $\delta$ with IC50 of 7 nM, 7 nM, 6 ...                                                                      |
| <b>Note</b>            | For research use only                                                                                                                                                                                     |
| <b>Formula</b>         | C <sub>26</sub> H <sub>26</sub> N <sub>4</sub> O <sub>3</sub>                                                                                                                                             |
| <b>MW</b>              | 442.5096                                                                                                                                                                                                  |
| <b>CAS Number</b>      | [133053-19-7]                                                                                                                                                                                             |
| <b>SMILES</b>          | <chem>O(C([H])([H])[H])C1C([H])=C([H])C2=C(C=1[H])C(C1C(N([H])C(C=1C1=C([H])N([H])C3=C([H])C([H])=C([H])C([H])=C13)=O)=O)=C([H])N2C([H])([H])C([H])([H])C([H])([H])N(C([H])([H])[H])C([H])([H])[H]</chem> |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                                                           |