

Product Datasheet

ML-228 (orb198443)

Description	ML-228 is a hIF pathway activator (EC50 values are 1.23 and 1.4 µM for HRE gene reporter assay and ...
Note	For research use only
Formula	C ₂₇ H ₂₁ N ₅
MW	415.49
CAS Number	[1357171-62-0]
SMILES	C1(C2=NC=CC=C2)=NC(NCC3=CC=C(C4=CC=CC=C4)C=C3)=C(C5=C C=CC=C5)N=N1
Expiration Date	12 months from date of receipt.