


Product Datasheet

BAF312 (Siponimod) (orb198453)

Description	BAF312 is a potent and selective agonist of S1P with EC50 value of 0.39nM for S1P1 receptors and 0.9...
Note	For research use only
Formula	C ₂₉ H ₃₅ F ₃ N ₂ O ₃
MW	516.59501
CAS Number	[1230487-00-9]
SMILES	N(CC1=CC=C(C(=N/OCC2=CC=C(C3CCCCC3)C(C(F)(F)F)=C2)/C)C=C1CC)1CC(C(O)=O)C1
Expiration Date	12 months from date of receipt.



Chemical structure of BAF312 (Siponimod)