


Product Datasheet

CD3 antibody (orb198485)

Description	Goat polyclonal antibody to CD3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CD3D (isoform A)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody has cross reactivity with Human
Protein Sequence	HETGRLSGAADTQ
Clonality	Polyclonal
MW	18.9; 14.5
Entrez	915
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat (Lane 1)...