


Product Datasheet

CD3 antibody (orb198490)

Description	Goat polyclonal antibody to CD3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CD3G
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody has cross reactivity with Human, Mouse, Rat, Dog, Pig, Cow
Protein Sequence	QDGVRSRASDKQT
Clonality	Polyclonal
MW	20.5
Entrez	917
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate U...