


Product Datasheet

NNMT antibody (orb198496)

Description	Goat polyclonal antibody to NNMT
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NNMT (aa171-182)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody has cross reactivity with Human
Protein Sequence	PDLPTYCRALRN
Clonality	Polyclonal
MW	29.6
Entrez	4837
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using NNMT antibody (orb198496).