


Product Datasheet

NNMT antibody (orb198497)

Description	Goat polyclonal antibody to NNMT
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NNMT
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody has cross reactivity with Human, Mouse, Rat
Protein Sequence	TSKDTYLSHFNP
Clonality	Polyclonal
MW	29.6
Entrez	4837
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 (Lane 1) a...