



Product Datasheet

APOF antibody (orb19907)

Description	Goat polyclonal antibody to APOF
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Apolipoprotein F
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 35kDa band observed in Human Liver lysates (calculated MW of 35.4kDa according to NP_001629.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	KSYDLDPGAGSLE
Clonality	Polyclonal
MW	35.4
Entrez	319
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 2.5-3.8µg/ml
Expiration Date	12 months from date of receipt.



Western blot anlysis of Human Liver Lys...



Immunohistochemical staining of Human Li...