

Product Datasheet

CHRNA3 antibody (orb19908)

Description	Goat polyclonal antibody to CHRNA3
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CHRNA3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx. 50kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 52.7kDa according to NP_000740.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	SRHVKKEHFISQ
Clonality	Polyclonal
MW	52.7
Entrez	1142
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19908
(1µg/ml)
staining of
Human
Cere...