

Product Datasheet

ZRANB1 antibody (orb19911)

Description	Goat polyclonal antibody to ZRANB1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TRABID / ZRANB1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx. 80kDa bands observed in Human Brain (Cerebellum) lysates (calculated MW of 81.0kDa according to NP_060050.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	KWLDRYRQIRPCTS
Clonality	Polyclonal
MW	81
Entrez	54764
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19911
(0.1ug/ml)
staining of
Human Ce...