


Product Datasheet

ABCD3 antibody (orb19913)

Description	Goat polyclonal antibody to ABCD3
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Target	ABCD3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. WB: Approx. 70kDa band observed in Human Kidney lysates (calculated MW of 75.5kDa according to NP_002849.1). Recommended concentration: 0.2-0.6 µg/ml.
Protein Sequence	PDGREDQKRKGISD
Clonality	Polyclonal
MW	75.5
Entrez	5825
Dilution Range	ELISA: 1:16000, WB: 0.2-0.6 µg/ml, IHC-P: 2.5-3.8 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...



Immunohistochemical staining of Human Te...