

Product Datasheet

OTUB2 antibody (orb19919)

Description	Goat polyclonal antibody to OTUB2
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	OTUB2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx. 27kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 27.2kDa according to NP_075601.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	YKTSHYNILYAADKH
Clonality	Polyclonal
MW	27.2
Entrez	78990
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19919
(1ug/ml)
staining of
Human
Cere...