

Product Datasheet

EGLN3 antibody (orb19921)

Description	Goat polyclonal to EGLN3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EGLN3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx. 30kDa band observed in Human Skeletal Muscle lysates (calculated MW of 27.7kDa according to NP_071356.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	RLGKYYVKERSK
Clonality	Polyclonal
MW	27.7
Entrez	112399
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19921
(0.1ug/ml)
staining of
Human Mu...