

Product Datasheet

TRPC4 antibody (orb19931)

Description	Goat polyclonal antibody to TRPC4
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TRPC4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 110kDa band observed in Rat Brain lysates (calculated MW of 112kDa according to Human NP_057263.1 and to Rat NP_536321.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	KHAKEEDSSIDYD
Clonality	Polyclonal
MW	112; 113; 103; 96.4; 94.9; 92.2
Entrez	7223
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 4.75 µg/ml
Expiration Date	12 months from date of receipt.



orb19931
(0.3ug/ml)
staining of
Rat Brai...