

Product Datasheet

CD22 Antibody (PE/Cy7) (orb1993554)

Description	CD22, also known as Siglec-2 (sialic acid-binding immunoglobulin-like lectin-2) is a transmembrane g...
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	human cell line Reh
Target	CD22
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P20273
Entrez	933
Expiration Date	12 months from date of receipt.