

Product Datasheet

CD11b Antibody (APC/Cy7) (orb1993593)

Description	CD11b (integrin alphaM subunit) is a 165-170 kDa type I transmembrane glycoprotein that non-covalent...
Reactivity	Human, Primate
Conjugation	APC/Cy7
Tested Applications	FC
Immunogen	Rheumatoid synovial cells and human monocytes.
Target	CD11b
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P11215
Entrez	3684
Expiration Date	12 months from date of receipt.