


Product Datasheet

APOB antibody (orb19937)


Description	Goat polyclonal antibody to APOB
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	APOB
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Bands corresponding to the expected size of ApoB100 were detected in Human serum samples from five different individuals. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	TDLHLRYQDKK
Clonality	Polyclonal
MW	516
Entrez	338
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 2 µg/ml
Expiration Date	12 months from date of receipt.



orb19937
(0.5ug/ml)
staining of
Human Se...



orb19937
(2µg/ml)
staining of
paraffin ...



orb19937
(2µg/ml)
staining of
paraffin ...