



Product Datasheet

CD207 Antibody (PE) (orb1993740)


Description	CD207 is a 40 kDa type II transmembrane glycoprotein of the C-type lectin family. It binds to mannos...
Reactivity	Human, Primate
Conjugation	PE
Tested Applications	FC
Immunogen	Fusion protein of human CD207 extracellular part and IgG Fc fragment
Target	CD207
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Extracellular and intracellular staining.
Isotype	Mouse IgG1 lambda
Clonality	Monoclonal
Uniprot ID	Q9UJ71
Entrez	50489
Expiration Date	12 months from date of receipt.



Flow cytometry multicolor surface staini...



Separation of human CD207 positive CD11c...



Flow cytometry surface staining pattern ...