

Product Datasheet

TCR Valpha24-Jalpha18 Antibody (PE) (orb1993742)

Description

TCR alpha V24-J18 chain is a marker of invariant NK

Reactivity

Human, Primate

Conjugation

PE

Tested

FC

Applications

Target

TCR Valpha24-Jalpha18

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype


Mouse IgG1 kappa

Clonality


Monoclonal

Expiration Date


12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation
of human
CD3
positive
TCR Val...



Flow
cytometry
surface
staining
pattern ...