



Product Datasheet

MGLL antibody (orb19939)


Description	Goat polyclonal antibody to MGLL
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	Monoglyceride Lipase
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.5-1.5 µg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 30kDa band observed in Human Adipose and Brain (Frontal Cortex) lysates (calculated MW of 33.3kDa according to NP_001003794.1). In transfected HEK293 transiently expressing MGLL a band of approx. 36kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	QDLPHLVNADGQY
Clonality	Polyclonal
MW	34.3; 33.3; 31.1
Entrez	11343
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Frontal Cortex lysates.



Immunohistochemical staining of Human Liver tissue.



Immunohistochemical staining of Human Liver tissue.