


Product Datasheet

TMPRSS4 antibody (orb19942)

Description	Goat polyclonal antibody to TMPRSS4
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TMPRSS4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 48kDa band observed in Human Bone Marrow lysates (calculated MW of 48.2kDa according to NP_063947.1 and of 47.6kDa according to NP_001077416.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	DSTRCNADDAYQGE
Clonality	Polyclonal
MW	48.2; 47.6
Entrez	56649
Dilution Range	ELISA: 1:32000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19942
(1ug/ml)
staining of
Human
Bone...