

Product Datasheet

ABHD12 antibody (orb19957)

Description	Goat polyclonal antibody to ABHD12
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ABHD12
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 40kDa band observed in Mouse Brain and Rat Brain lysates (calculated MW of 45.1kDa according to Human NP_001035937.1 and of 45.3kDa according to Mouse NP_077785.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	REFLGKSEPEHQH
Clonality	Polyclonal
MW	45.1; 45.6
Entrez	26090
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19957
(0.3ug/ ml)
staining of
Mouse B...