


Product Datasheet

AADAT antibody (orb19961)

Description	Goat polyclonal to AADAT
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	AADAT
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 45kDa band observed in Human Liver lysates (calculated MW of 47.4kDa according to NP_057312.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	KPEDAKNPQKNTPK
Clonality	Polyclonal
MW	47.4
Entrez	51166
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19961
(1ug/ml)
staining of
Human Live...