


Product Datasheet

SPI1 antibody (orb19964)


Description	Goat polyclonal antibody to SPI1
Species/Host	Goat
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PU.1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 40kDa band observed in lysates of DAUDI and K562 (calculated MW of 31.2kDa according to NP_001074016.1 and of 31.1kDa according to NP_003111.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Harendza S et al, J Biol Chem. 2000 Jun 30;275(26):19552-9.; PMID: 10867017). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	DLYQRQTHEYY
Clonality	Polyclonal
MW	31.2; 31.1
Entrez	6688
Dilution Range	ELISA: 1:32000, WB: 0.03-0.1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb19964 (0.03ug/ml)
staining of DAUDI l...



Immunohistochemical
staining of Human
Co...



Immunohistochemical
staining of Human
Sp...