

Product Datasheet

FGFR2 antibody (orb19966)

Description	Goat polyclonal antibody to FGFR2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FGFR2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 100kDa band observed in Human lysates of cell line A549 (calculated MW of 92.0kDa according to NP_000132.1 and of 92.1kDa according to NP_075259.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	GREKEITASPDY
Clonality	Polyclonal
MW	92.0; 92.1; 86.2; 79.8; 79.3; 79.4; 79.2; 76.4
Entrez	2263
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19966
(0.3ug/ml)
staining of
A549 lys...