

Product Datasheet

MYL3 (myosin, light chain 3, alkali; ventricular, skeletal, slow) Blocking Peptide (the N terminal of protein) (orb1997439)

Description	MYL3 (myosin, light chain 3, alkali; ventricular, skeletal, slow) Blocking Peptide (the N terminal o...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Uniprot ID	Q5R887
NCBI	NP_000249
Expiration Date	6 months from date of receipt.