

Product Datasheet

DMP1 antibody (orb19977)

Description	Goat polyclonal to DMP1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	DMP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx. 55kDa band observed in Human Ileum and Kidney lysates (calculated MW of 55.8kDa according to NP_004398.1 and of 54.1 according to NP_001073380.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	ENSEQANTQDSG
Clonality	Polyclonal
MW	55.8; 54.1
Entrez	1758
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19977
(1µg/ml)
staining of
Human Kidn...