


Product Datasheet

IKZF1 antibody (orb19989)

Description	Goat polyclonal antibody to IKZF1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IKZF1 / IKAROS
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 65kDa band observed in lysates of cell line MOLT4 (calculated MW of 57.5kDa according to NP_006051.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	DLLRAASENSQDALR
Clonality	Polyclonal
MW	57.5; 52.7; 48.3; 47.2; 47.6; 43.5; 42.4; 41.2
Entrez	10320
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of MOLT4 lysate us...