

Product Datasheet

MUPCDH antibody (orb20000)

Description

Goat polyclonal to MUPCDH

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Target

MUPCDH

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

This antibody is expected to recognise both reported isoforms (NP_068743.2; NP_112554.2).ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Preliminary experiments in Human Colon, Ileum and Kidney lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates

Protein Sequence

PYDAPGGDDSY

Clonality

Polyclonal

MW

88.2; 69.6

Entrez
53841
Dilution Range

ELISA: 1:8000, WB: 1 µg/ml

Expiration Date

12 months from date of receipt.