



Product Datasheet

Apom antibody (orb20015)

Description	Goat polyclonal antibody to Apom
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Apolipoprotein M (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.
Protein Sequence	KGNMELRTEGRPD
Clonality	Polyclonal
MW	21.3
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb20015 (0.5ug/ml)
staining of Human
PL...



Immunohistochemical
staining of Human
Li...