



Product Datasheet

GPR120 antibody (orb20019)

Description	Goat polyclonal antibody to GPR120
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	GPR120
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Preliminary experiments in Human Colon and Duodenum lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	YSESHQIRVSQQD
Clonality	Polyclonal
MW	42.2
Entrez	338557
Dilution Range	ELISA: 1:16000, WB: 1µg/ml, IHC-P: 2-4 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Co...



Immunohistochemical staining of Human Co...