

Product Datasheet

LDLRAD2 Peptide - middle region (orb2002402)

Description	LDLRAD2 Peptide - middle region
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	26kDa
Uniprot ID	Q5SZI1
NCBI	NP_001013715
Expiration Date	6 months from date of receipt.