


Product Datasheet

LTF antibody (orb20033)

Description	Goat polyclonal to LTF
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	LTF
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded blood vessel shows strong staining of cytoplasm of neutrophils.. Recommended concentration, 2-4ug/ml.WB: Approx. 85kDa band observed in Human Peripheral Mononucleocyte lysates (calculated MW of 78.2kDa according to NP_002334.2). Recommended concentration: 0.05-0.2 µg/ml.
Protein Sequence	ENYKSQQSSDPDP
Clonality	Polyclonal
MW	78.2
Entrez	4057
Dilution Range	ELISA: 1:32000, WB: 0.05-0.2 µg/ml, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Peripheral Bloo...



immunohistochemical staining of paraffin...