

Product Datasheet

COG7 antibody (orb20034)

Description	Goat polyclonal antibody to COG7
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	COG7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 85kDa band observed in lysates cell lines HeLa, Jurkat and MOLT4 (calculated MW of 86.3kDa according to NP_705831.1). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	KTRPEDYRQVSK
Clonality	Polyclonal
MW	86.3
Entrez	91949
Dilution Range	ELISA: 1:32000, WB: 0.03-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb20034
(0.03ug/ml)
staining of
MOLT4 I...