


Product Datasheet

Fto antibody (orb20045)

Description	Goat polyclonal antibody to Fto
Species/Host	Goat
Reactivity	Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FTO (Mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx. 60kDa band observed in Human Pancreas lysates (calculated MW of 58.0kDa according to NP_036066.2). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	QQKPDCRPYWEKDD
Clonality	Polyclonal
MW	58
Entrez	79068
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



orb20045
(0.03ug/ml)
staining of
Human P...