


Product Datasheet

APOL1 antibody (orb20048)

Description	Goat polyclonal antibody to APOL1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	APOL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx. 45kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 44.0kDa according to NP_003652.2;). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NNNYKILQADQE
Clonality	Polyclonal
MW	44.0; 45.9
Entrez	8542
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Placenta ...