

Product Datasheet

41889 antibody (orb20054)

Description	Goat polyclonal antibody to 41889
Species/Host	Goat
Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SEPT7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 45kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 50.7kDa according to NP_001779.3 and of 55.6kDa according to NP_001011553.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	DNNKNKGQLTKSP
Clonality	Polyclonal
MW	50.7; 50.6; 46.6
Entrez	989
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20054
(0.5ug/ml)
staining of
nuclear ...