


Product Datasheet

KRT13 antibody (orb20058)

Description	Goat polyclonal antibody to KRT13
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	KRT13
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.IHC: In paraffin embedded Human Tonsil shows strong layered, filamentous staining. Recommended concentration, 2-4ug/ml.WB: Approx. 55kDa band observed in Human Lung, Skeletal Muscle and Skin lysates (calculated MW of 49.6kDa according to NP_705694.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Sajjan, Microbiology. 2003 Apr;149(Pt 4):961-71.; PMID: 12686638). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	TSNASGRRTSDVRRP
Clonality	Polyclonal
MW	49.6; 45.8
Entrez	3860
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



orb20058
(0.1ug/ml)
staining of
Human Lu...



orb20058
(2.5ug/ml)
staining of
paraffin...