

Product Datasheet

PLEKHB1 antibody (orb20063)

Description	Goat polyclonal antibody to PLEKHB1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PHR1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx. 25kDa band observed in Rat Brain lysates (calculated MW of 25.2kDa according to Rat NP_742030.1 (Evectin-1) and of 27.2kDa according to Human NP_067023.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	PPDSALESPFEE
Clonality	Polyclonal
MW	27.2; 23.2
Entrez	58473
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20063
(1ug/ml)
staining of
Rat Brain ...