


Product Datasheet

APOL2 antibody (orb20086)

Description	Goat polyclonal antibody to APOL2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	APOL2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx. 37.1kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 38kDa according to NP_112092.1; NP_663612.1). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	TKIHEMLQPGQD
Clonality	Polyclonal
MW	37.1
Entrez	23780
Dilution Range	ELISA: 1:64000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



orb20086
(0.03ug/ml)
staining of
Human F...